FOSTER[®] PATIENT RECLINER

Information Guide



ARMS

Pivoting armrests on both sides of the chair allow a side access for the patient and caregiver. Arm caps available in grey polyurethane to match other controls.

COMFORT

Passive built-in lumbar support, dual density seat suspension with molded foam, sinuous springs and an optional adjustable headrest all speak to the user needs and provide comfort.



CLEANABILITY

The Foster Patient Recliner offers a 360° wipeout along with standard removable covers and a removable seat.

BACK RECLINE AND FOOTREST

The back and footrest work independently, with the ability to lay completely flat, which offers a comfortable sleeping alternative.

CONTROLS

The Patient Recliner can be configured with various controls on each side, including a footrest lever, central brake and steer, recline activation levers and reverse incline.